

## KIMTECH Science™ Kimwipes™ Delicate Task Wipers

### Product Description:

KIMTECH Science™ Kimwipes™ are delicate, light-duty wipers to handle precise tasks such as cleaning surfaces, parts, instruments in labs, and laboratory lenses. These wipers easily clean liquids, dust and small particles. Plus, the anti-static dispensing design reduces electrostatic discharge. These wipers are an efficient choice, thanks to the waste-reducing, one-at-a-time dispensing.

### Product Details:

Product Code	Color	Sheet Size	Number of Ply	Case Count
34155	White	215 x 110 mm	1	280 sheets/box 60 boxes/case
34120	White	215 x 110 mm	1	280 sheets/box 30 boxes/case

### Physical Properties:

The following targets are current as of the date of this revision, subject to change at any time without prior notice.

Physical property data is as produced; subsequent storage conditions and the age of the product, which are beyond Kimberly-Clark's control, may affect the properties of the product as received.

Characteristics	Target	Method
Basis Weight	18.7 g/m <sup>2</sup> 1ply	STM-00109
MDT	2200 g/76mm sheet	STM-00140
CDT	1450 g/76mm sheet	STM-00140
CD-WT	380 g/76mm sheet	STM-00140
Gelbo Lint	20 (>50 microns)	STM-00135
Silicone	0 ppm	Boeing IR

\*All care has been taken to ensure that the information contained herein is as accurate as possible at the time of publication, however may change in the lifetime of this material. Accordingly, the specifications for the products may be subjected to change.

### Workmanship Standard:

Wipers are made in accordance with reasonable industry practice with respect to holes, tears, wrinkles, cleanliness and foreign materials or dirt. Edges are to be cleanly cut and not ragged. Wipers are folded as specified. Wipers have no disagreeable odour, either wet or dry, in accordance with reasonable industry practice.

For additional information please contact your local Kimberly-Clark representative  
[www.kcprofessional.com](http://www.kcprofessional.com)